

HALUX N30



medicinelaser^{*}
LATAM 



Small but mighty! HALUX N30

The HALUX N30 is a solid examination fixture with uncompromisingly good light and a very attractive price. The illumination field has a diameter of 18 cm and an illuminance of up to 30,000 lux at 0.5 m which makes it perfectly equipped for the versatile daily routine in the clinic. The good colour rendering properties and the stepless dimming function offer high functionality.



HALUX N30	technical data
luminous power *	30,000 Lux @ 0.5 m
colour temperature *	4400 K
lighting field @ 1.0 m *	Ø 18 cm
colour rendering index	Ra > 93
usage	touch panel on luminaire head
class of protection	II
protection rating	IP 20

special features	
HALUX N30	<ul style="list-style-type: none"> continuously dimmable maintenance-free spring arm
variations	<ul style="list-style-type: none"> gynecology ambulance
fastening	<ul style="list-style-type: none"> pin universal roller stand
accessories	<ul style="list-style-type: none"> wall mount roller stand fasting hardware universal mount table clamp

HALUX N30 – variations



Variety of versions

The HALUX N30 is available in different versions. Depending on your requirements, it can be mounted directly on the wall or fixed to a rail mount with a rail clamp. The ambulance version can be mounted on C-rails. The gynecology version is suitable for connection to the gynecology chair.

Conceivably simple

The 18 cm large illumination field can be easily controlled via the touch panel on the luminaire head. The illuminance can be dimmed steplessly to adapt to the lighting requirements. Cleaning is also simple. This is ensured by the seamless design.

HALUX N30 P SGH and SGV

pin version | flex arm with foot joint



HALUX N30 P SH and SV

pin version | flex arm



HALUX N30 – variations



HALUX N30 P F1

pin version | spring arm with foot joint



HALUX N30 P Amb

pin version | head joint | ambulance



HALUX N30 F

integrated roller stand



HALUX N50 and N30 – gynecology





Lighting solutions for gynecology & urology

We support gynecological and urological clinics as well as proctologists with optimal light. Very good colour rendering properties with Ra values > 93 ensure a natural representation of visible pathologies, tissue and vessels. Unobtrusive spotlights on the examination chair provide ideal light for patient examination.

HALUX N50-3 P SX Gyn

pin version | flex arm | gynecology

- 50,000 Lux @ 0.5 m
- optional: 5-stage dimmable
- adjustable colour temperature (3300/ 3800/ 4400 K)



Securely fastened
rail clamp with brake function

HALUX N30-1 P SGV Gyn

pin version | flex arm | gynecology

- 30,000 Lux @ 0.5m
- continuously dimmable illumination field with head connection vertical/horizontal
- colour temperature fixed 4400 K



Optimal vision
thanks to homogeneous, spot-like illumination field

medicinelaser
LATAM

www.medicinelaser.com



ARGENTINA
Av. Doctor Honorio
Pueyrredón 707, Caballito,
CABA, Argentina.
M: +54 9 11 3404-7158

CHILE
Sucre 1355, Ñuñoa,
Santiago, Chile.
O: +56 2 2269 1981
M: +56 9 4409 7105

ECUADOR
Av. República 1058,
Edificio La Fontana,
Quito, Ecuador.
M: +593 99 120 4597

PERÚ
Calle Almirante Lord
Cochrane 322,
Miraflores, Lima, Perú.
M: +51 983 552 813

+ OFICINAS CENTRALES

ventas@medicinelaser.com