

HALUX N50



Versatile classic. HALUX N50

HALUX N50 is the right luminaire for anyone looking for good lighting quality. The illuminance of up to 50,000 lux at 0.5 m and the 18 cm large light field provide optimal examination conditions. It is versatile for examinations in dermatology, intensive care, neonatology, gynecology, obstetrics or urology, but also at the patient's bedside or for the preparation of a surgical procedure.



HALUX N50	technical data
luminous power *	50,000 Lux @ 0.5 m
colour temperature *	fixed: 4400 K adjustable: 3300 K / 3800 K / 4400 K **
lighting field @ 0.5 m *	Ø 18 cm
colour rendering index	Ra > 93
usage	touch panel on luminaire head or rocker switch (I/O)
class of protection	II
protection rating	IP 20 (luminaire head)

special features	
HALUX N50	<ul style="list-style-type: none"> • 3 colour temperatures ** • 5-stage dimmable ** • maintenance-free
variations	gynecology
fastening	pin universal
accessories	<ul style="list-style-type: none"> • wall mount • roller stand • fastening hardware • universal mount • table clamp

HALUX N50 – variations



HALUX N50 P FX

pin version | spring arm



HALUX N50 P SX

pin version | flexarm



Large assortment

HALUX N50 impresses with a well thought out variety: with a robust flexible arm or a positionally stable spring arm, LEDs with defined or adjustable colour temperatures and various mounting options. A special version for connection to the examination chair in gynecology and proctology is also available.

Durable use

High-quality materials and long-life LEDs ensure a flawless use. The closed design facilitates cleaning and prevents the penetration of cleaning agents into the luminaire head.

Product highlights



Intuitive operation
simple and within easy reach via touch panel



Individual light intensity
thanks to 5-stage dimming the light is easily adjustable



Perfect colour rendering
optimum light colour for precise and accurate diagnoses



Optimal contrast vision
optional: 3 colour levels for detail fidelity as well as perfect contrasts

3300 K | 3800 K | 4400 K



medicinelaser
L A T A M

www.medicinelaser.com



ARGENTINA
Av. Doctor Honorio
Pueyrredón 707, Caballito,
CABA, Argentina.
M: +54 9 11 3404-7158

CHILE
Sucre 1355, Ñuñoa,
Santiago, Chile.
O: +56 2 2269 1981
M: +56 9 4409 7105

ECUADOR
Av. República 1058,
Edificio La Fontana,
Quito, Ecuador.
M: +593 99 120 4597

PERÚ
Calle Almirante Lord
Cochrane 322,
Miraflores, Lima, Perú.
M: +51 983 552 813

+ OFICINAS CENTRALES

ventas@medicinelaser.com