

OPTICLUX Hand



medicinelaser
LATAM 



Strikingly different. OPTICLUX Hand

The OPTICLUX Hand magnifier luminaire with ergonomic handle provides support for demanding visual tasks. This is ensured by the large dimmable light field in daylight quality, excellent colour rendering properties and long battery life. Dermatologists and veterinarians can look forward to the powerful Wood's Light, which is available as an option in combination with the daylight mode.



OPTICLUX Hand	technical data
luminous power *	8,000 Lux @ 0.15 m
colour temperature *	<ul style="list-style-type: none"> 6500 K optional: 365 nm (Wood's Light)
lens dimension	Ø 16 cm
colour rendering index	Ra > 95
usage	touch panel on handle
class of protection	II
protection rating	IP 20



special features
<ul style="list-style-type: none"> dimmable rechargeable battery with USB-C charging interface
OPTICLUX Hand <ul style="list-style-type: none"> Wood's Light mode (optional) lens cover 1.9x magnification
variations <ul style="list-style-type: none"> battery operated handheld version

OPTICLUX Hand – variations



OPTICLUX Hand

Daylight & Wood's Light



OPTICLUX Hand + Table charger

Setprice available



First class light

The OPTICLUX hand magnifier light delivers unique light output of 8,000 lux at a distance of 15 cm. Excellent colour rendering properties > Ra 95 ensure ideal vision.

Large field of view

The scratch-resistant plastic lens with 160 mm diameter and 3.5 diopters is ideal for full-body examinations and for assessing dermatological diseases on patients.

Product highlights



Optional with Wood's Light

two lighting modes in one device: optional Wood's Light for examination of fungal or bacterial infection on the skin or scalp and pigmentary changes

For detailed description of Wood's Light see page 110.



Ergonomic

optimal handling for right- and left-handed



Table charging station
for storage and charging, white powder-coated with a sturdy aluminium base



Easy to operate
touch panel on the handle with comfortable, stepless dimming



Battery operated
up to 3 hours battery runtime, flexible charging through USB-C interface

medicinelaser
LATAM

www.medicinelaser.com



ARGENTINA
Av. Doctor Honorio
Pueyrredón 707, Caballito,
CABA, Argentina.
M: +54 9 11 3404-7158

CHILE
Sucre 1355, Ñuñoa,
Santiago, Chile.
O: +56 2 2269 1981
M: +56 9 4409 7105

ECUADOR
Av. República 1058,
Edificio La Fontana,
Quito, Ecuador.
M: +593 99 120 4597

PERÚ
Calle Almirante Lord
Cochrane 322,
Miraflores, Lima, Perú.
M: +51 983 552 813

+ OFICINAS CENTRALES

ventas@medicinelaser.com