

TRIANGO 60

The compact one that has what it takes



TRIANGO 60 – your new treatment & minor surgery light



TRIANGO 60 C
ceiling version




TRIANGO 60 W
wall mount



TRIANGO 60 F
roller stand



TRIANGO 60 F Battery
roller stand



TRIANGO 60	technical data
Fitted with	1x LED-Modul
Luminous power *	60 000 Lux / 1,0 m
Colour temperature *	fixed value 4 700 K
TRIANGO 60-1 F	changeable value
TRIANGO 60-3 F	3 700 / 4 300 / 4 700 K
Lighting field 1,0 m*	ø 18 cm
Depth of illumination *	140 cm
Colour rendering index *	Ra > 94
Operation equipment	Electronic converter
Voltage	100-240 V; 50/60 Hz
Luminaire body	Polyamide (PA)
Power consumption	< 30 W
Usage	Touch panel on cardanic system
Technology	Single-pole secondary
Balance of articulation joints	Spring (weight compensation)
Class of protection	I
Protection rating	IP 43 (illumination head) IP 20 (UPSmedv Pack)
Standards met	EN 60601-1, EN 60601-2-41

Special features	
TRIANGO 60-1	<ul style="list-style-type: none"> Continuously dimmable
TRIANGO 60-3	<ul style="list-style-type: none"> Continuously dimmable 3 colour temperatures
TRIANGO Fokus 60-1	<ul style="list-style-type: none"> Continuously dimmable Focusable
TRIANGO 60-1 F Battery	<ul style="list-style-type: none"> Continuously dimmable UPSmed Pack
Mounting	<ul style="list-style-type: none"> Ceiling mount Wall mount Roller stand
Cord length	
TRIANGO 60-1 C	Direct power supply
TRIANGO 60-1 W	Direct power supply
TRIANGO 60-1 F	3 m
TRIANGO 60-1 F Battery	3 m
Weight (net)	
TRIANGO 60-1 C	13 kg
TRIANGO 60-1 W	12,5 kg
TRIANGO 60-1 F	19 kg
TRIANGO 60-1 F Battery	31 kg
Accessories	<ul style="list-style-type: none"> Sterilizable handle Disposable cover for handle

Versatile
versions with focusable or colour changeable illumination field

Homogeneous light field
thanks to unique diamond optic

Easy cleaning
streamlined design and the removable handle

Hygiene made simple
good design does without unnecessary details and attaches importance to closed building structures

Can be used everywhere
TRIANGO 60 is suitable for the versatile daily routine in family practices and clinics as well as in gynecology, dermatology or veterinary medicine



Ergonomic handling
easy alignment for right- and left-handed



Intuitive operation
thanks to an intuitive touch panel on the cardanic, the illumination from 60 000 lux at 1 m can be infinitely adjusted

Order number: TRIANGO 60-1	
TRIANGO 60-1 C	D16.250.000-00813741 D16.253.000-00814312 - USA
TRIANGO Fokus 60-1 C	D16.252.000-00814311 D16.255.000-00814314 - USA
TRIANGO 60-1 W	D16.256.000-00813742 D16.274.000-00814327 - USA
TRIANGO Fokus 60-1 W	D16.258.000-00814316 D16.276.000-00814329 - USA
TRIANGO 60-1 F	D16.262.000-00813744 - CH D16.259.000-00813743 - Schuko D16.271.000-00813747 - UK D16.265.000-00813745 - USA
TRIANGO Fokus 60-1 F	D16.264.000-00814320 - CH D16.261.000-00814318 - Schuko D16.273.000-00814326 - UK D16.270.000-00814322 - USA
TRIANGO 60-1 F Battery	D16.279.000-00814482 - CH D16.278.000-00814481 - Schuko D16.282.000-00814485 - UK D16.280.000-00814483 - USA

Order number: TRIANGO 60-3	
TRIANGO 60-3 C	D16.251.000-00814310 D16.254.000-00814313 USA
TRIANGO 60-3 W	D16.257.000-00814315 D16.275.000-00814328 USA
TRIANGO 60-3 F	D16.263.000-00814319 - CH D16.260.000-00814317 - Schuko D16.272.000-00814325 - UK D16.266.000-00814321 - USA



* -10% / +20% tolerance



Headquarters

Derungs Licht AG
Hofmattstrasse 12
9200 Gossau
Switzerland
+41 71 388 11 66
mailbox@derungs.swiss
www.derungs.swiss

medicinelaser
L A T A M

www.medicinelaser.com



ARGENTINA
Av. Doctor Honorio
Pueyrredón 707, Caballito,
CABA, Argentina.
M: +54 9 11 3404-7158

CHILE
Sucre 1355, Ñuñoa,
Santiago, Chile.
O: +56 2 2269 1981
M: +56 9 4409 7105

ECUADOR
Av. República 1058,
Edificio La Fontana,
Quito, Ecuador.
M: +593 99 120 4597

PERÚ
Calle Almirante Lord
Cochrane 322,
Miraflores, Lima, Perú.
M: +51 983 552 813

+ OFICINAS CENTRALES

ventas@medicinelaser.com



Derungs Licht AG is
a company of the
Waldmann Group

www.derungs.swiss



Register and benefit:
+ 1 year warranty and
exclusive offers

www.derungs.swiss/product-registration