

# TRIANGO 80



medicinelaser<sup>\*</sup>  
LATAM 



reddot winner 2012

## A proven classic. TRIANGO 80

With 80,000 lux at a distance of one meter and a light field diameter of 18 cm, the medical staff has a powerful treatment light at its disposal.

The illumination field is pleasantly homogeneous and continuously dimmable.

The arm system enables simple and accurate alignment of the lamp head position.

The handle makes it easy for right-handed and left-handed users to position the light.



TRIANGO 80	technical data
luminous power *	80,000 Lux @ 1.0 m
colour temperature *	ca. 4500 K
lighting field @ 1.0 m *	Ø 18 cm
depth of illumination	124 cm
colour rendering index	Ra > 95
usage	rotary switch on the luminaire head
class of protection	I
protection rating	IP 43 (luminaire head) IP 20 (UPSmed pack)

special features	
TRIANGO 80-1	<ul style="list-style-type: none"> <li>continuously dimmable</li> </ul>
fastening	<ul style="list-style-type: none"> <li>ceiling</li> <li>wall</li> <li>roller stand</li> <li>roller stand UPSmed pack</li> </ul>

# TRIANGO 80 – variations





TRIANGO 80 C ceiling version



TRIANGO 80 W wall mount



TRIANGO 80 F roller stand



TRIANGO 80 F Battery-UPSmed Pack



medicinelaser  
LATAM

[www.medicinelaser.com](http://www.medicinelaser.com)



ARGENTINA  
Av. Doctor Honorio  
Pueyrredón 707, Caballito,  
CABA, Argentina.  
M: +54 9 11 3404-7158

CHILE  
Sucre 1355, Ñuñoa,  
Santiago, Chile.  
O: +56 2 2269 1981  
M: +56 9 4409 7105

ECUADOR  
Av. República 1058,  
Edificio La Fontana,  
Quito, Ecuador.  
M: +593 99 120 4597

PERÚ  
Calle Almirante Lord  
Cochrane 322,  
Miraflores, Lima, Perú.  
M: +51 983 552 813

+ OFICINAS CENTRALES

[ventas@medicinelaser.com](mailto:ventas@medicinelaser.com)